


Easy-load Pump Head



YZ1515x, YZ2515x

Classic design, simple operation, wide flow range.



PPS housing
(Excellent chemical resistance)





YZII15, YZII25

Patented Design

Patent Number: 200620026529.1

Based on the classic design of YZ1515X/YZ2515X, YZII15/YZII25 have a special tubing retention linkage which makes tubing loading easily and rapidly. YZII25 accepts more suitable tubing than YZ2515X and provides a wider flow range.



Tubing Loading Comparison



Tubing clamber needs to be lifted manually
YZ1515X/YZ2515X



Tubing clamber can be automatically opened
YZII15/YZII25

Model	Product Code	Roller Number	Roller Material	Housing Material	Speed (rpm)	Tubing	Max. Flow Rate (mL/min)	Weight (kg)
YZ1515x-A	05.01.52A	3	Stainless steel	PESU	< 600	13° 14° 19° 16° 25° 17° 18°	2200	0.4
YZ1515x-B	05.01.52B	6						
YZ1515x-C	05.01.52C	3						
YZ1515x-D	05.01.52D	6						
YZ2515x-A	05.01.53A	3		PESU		15° 24°		
YZ2515x-C	05.01.53C		PPS					
YZII15	05.01.55A	3	Stainless steel	PESU	< 600	13° 14° 19° 16° 25° 17° 18°	2200	0.35
YZII25	05.01.56A							